

# Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **Regenerier-Tabts**

Page 1/5

Revised at: 08.11.2007

Print date: 13.11.2007

## 1. Identification of the substance/preparation and of the company/undertaking

Trade name: **Regenerier-Tabts**

Recommended use: Special tablets for cleaning and descaling of the Siemens SIRONA® Hygiene Centers and the SIROCLAVES / SIROMINI.

Effect of the preparation: Cleansing and descaling agent

Distributed by: Sirona Dental Systems GmbH  
Fabrikstr. 31  
D-64625 Bensheim  
Telephone: +49 6251 161616

Manufacturer: ALPRO MEDICAL GMBH  
Mooswiesenstr. 9  
D-78112 St. Georgen (Germany)  
Telephone: +49 7725 9392-0  
Telefax: +49 7725 9392-91  
Email: alpro@alpro-medical.de  
Internet: www.alpro-medical.com

Technical contact point: +49 7725 9392-0

Emergency phone: +49 7725 9392-0 (business hours)  
+49 761 19240 Poisoning-information-central, Freiburg (24h emergency call)

## 2. Hazards identification

R36/38 Irritating to eyes and skin.

## 3. Composition/information on ingredients

Chemical characterization: Mixture of substances listed below with nonhazardous additions.

<i>Hazardous ingredients</i>	<i>CAS-No.</i>	<i>Weight-%</i>	<i>Code letter</i>	<i>R-phrases*</i>
Citric acid anhydrate	77-92-9	15-30%	Xi	36
Sodium-cumolsulfonate	28348-53-0	15-30%	Xi	36/38
Oleoylmethyltaurine-sodium-salt	137-20-2	5-15%	Xi	41
Amidosulfonic acid	5329-14-6	5-15%	Xi	36/38; 52/53
Trisodiumnitrilotriacetate	5064-31-3	1-5%	Xn	22, 36

\*Full text of R-phrases: see under section 16

## 4. First aid measures

General information: Note precautions for self-protection.

On inhalation: Move affected person into fresh air and keep still and warm. In case of irritation of the airways, seek medical advice.

On contact with skin: Wash off immediately with soap and water. In case of skin reactions, seek medical advice.

On contact with the eyes: Rinse the eyes with water with the eyelids open for a sufficient length of time, then seek ophthalmological advice immediately.

If swallowed: Rinse out mouth and then drink plenty of water. Do not induce vomiting. Seek medical advice.

# Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **Regenerier-Tabts**

Page 2/5

Revised at: 08.11.2007

Print date: 13.11.2007

Information for the doctor: First aid, decontamination, symptomatic treatment.

---

## 5. Fire-fighting measures

Suitable extinguishing media:	Water, foam, dry fire-extinguishing media, carbon dioxide
Extinguishing media which must not be used for safety reasons:	None known
Special exposure hazards arising from the substance or preparation itself, combustion products, resulting gases:	None known
Special protective equipment for firefighting:	Not necessary
Additional information:	Cool endangered containers with water spray jet and if possible remove off fire area.

---

## 6. Accidental release measures

Personal precautions:	Use personal protective clothing.
Environmental precautions:	Pick up dry.
Methods for cleaning up/taking up:	Treat recovered material as prescribed in: Disposal considerations.

---

## 7. Handling and storage

Hints for safe handling:	Avoid contact with skin and eyes.
Hints for protection against fire and explosion:	Product is not flammable.
Further informations:	Not necessary
Requirements for storage rooms and vessels:	Keep tightly closed and store only in original container.
Hints on storage assembly:	Not necessary
Further informations on storage conditions:	Store only in original container at temperatures not exceeding 60 °C. Protect from humidity.
Storage class:	LGK 11

---

## 8. Exposure controls / Personal protection

### Constituents with occupational exposure limits resp. biological occupational exposure limits requiring monitoring:

Does not contain substances above concentration limits fixing an occupational limit.

### General health and safety measures:

- Do not eat, drink, smoke while working.
- Wash hands before breaks and on finishing work.
- Preventive skin protection using barrier ointment.
- Keep away of foods and drinks.
- Avoid contact with eyes, skin and clothing.
- Avoid inhalation of aerosols and vapours.

# Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **Regenerier-Tabts**

Page 3/5

Revised at: 08.11.2007

Print date: 13.11.2007

## Personal protection equipment:

Respiratory protection:	Not necessary if handled cautiously.
Hand protection:	E.g. natural rubber- or PVC-gloves.
Eye protection:	Tightly sealed goggles.
Protective clothing:	Not necessary if handled cautiously.

## 9. Physical and chemical properties

Form:	tablet	
Colour:	white-slightly yellowish	
Odour:	acidic	
pH-value (2 g/l H <sub>2</sub> O):	3,0-4,0	(20°C)
Boiling point/Boiling range:	not relevant to safety	
Melting point/Melting range:	not relevant to safety	
Flash point (DIN, TRbF):	>55°C	
Flammability (solid, gases):	not determinable	
Explosive properties:	none	
Lower explosion limit:	none	
Upper explosion limit:	none	
Ignition temperature:	not necessary	
Oxidizing properties:	none	
Vapour pressure:	not necessary	at ...°C
Density:	none	at ...°C
Bulk density:	not determinable	at ...°C
Solubility in water:	soluble	
in ethanol:	soluble	
in hexan:	not soluble	
Partition coefficient:	not necessary	(n-octanol/water)
Viscosity (sort):	not necessary	at ...°C
Electrical conductivity (2 g/l H <sub>2</sub> O):	1000-1500 µS/cm	(20°C)
Refractive index n <sub>D</sub> :	none	
Optical rotation α <sub>D</sub> :	not necessary	

## 10. Stability and reactivity

Conditions to avoid:	No hazardous decomposition products if stored and handled as prescribed.
Materials to avoid:	No hazardous reactions with other substances known.
Hazardous decomposition products:	No decomposition if used as prescribed.

# Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **Regenerier-Tabts**

Page 4/5

Revised at: 08.11.2007

Print date: 13.11.2007

## 11. Toxicological information

Acute toxicity:	No acute toxicity known.
Primary irritant effect on the skin:	Irritating to skin.
Primary irritant effect on the eye:	Irritating to eyes.
Mutagenic potential:	No mutagenic potential known.
Sensitization:	No sensitizing effects known.

The product was classified on the basis of the calculation procedure of the preparation directive 1999/45/EC.

Based on manufacturer data of the main components the acute toxicity, the skin irritation, the mucous membrane irritation and the mutagenic potential of the preparation were evaluated. However, some data are not complete regarding particular main components. Nevertheless according to the experience of the manufacturer there are no other hazards expected than those which are already mentioned on the label.

## 12. Ecological information

Ecotoxicity:	Inhibition of degradation activity in activated sludge is not to be anticipated during correct introduction of low concentrations.
Persistence and degradation:	Ingredients are biodegradable.
Other adverse effects:	No other adverse effects known.
Further information:	Do not allow product to reach ground water or waters and do not empty undiluted product into sewers.

## 13. Disposal considerations

### Disposal of the product:

In the EC member states there are no standardized regulations for the disposal of chemicals. In Germany the law of the cycle of economy and waste (KrW/AbfG) prescribes the principle of utilization and therefore a distinction must be made between "waste for utilization" and "waste for removal". Furthermore special cases – particularly in case of delivery – are also regulated by the federal states. Please contact the corresponding office (authority or waste disposal company) who will inform you about the utilization or removal.

Waste code no.\*: not specified

### Disposal of uncleaned as well as completely emptied packings:

Disposal according to official regulations. Contaminated packings are to be treated like the substance.

Waste code no.\*: not specified

### Further information:

If there is no official regulation, non-contaminated packings can be treated like domestic waste or can be recycled.

\*Waste code no.: Waste incurred is classified following the code of the European waste register according to kind of waste and branch.

# Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **Regenerier-Tabs**

Page 5/5

Revised at: 08.11.2007

Print date: 13.11.2007

## 14. Transport information

No hazardous material as defined by the transport regulations (ADR/RID, IMDG-Code, ICAO-TI/IATA-DGR)

## 15. Regulatory information

### Labelling

Code letter:	Xi	
Hazard designation:	Irritant	
contains:	-	
R-phrases:	R36/38	Irritating to eyes and skin.
S-phrases:	S26	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
	S37/39	Wear suitable gloves and eye/face protection.

Special labelling of certain preparations (preparation directive 1999/45/EC annex V):  
not necessary

### National regulations (Germany)

Restrictions of occupation:	Juvenile in relation to § 22JArbSchG
Störfallverordnung:	-
Classification according to VbF:	Not subject to the VbF
Technische Anleitung Luft:	-
Water hazard class:	WGK 2 (water hazardous); self-classification

## 16. Other information

### Full text of R-phrases appearing in section 2:

(This is not the classification of the product Regenerier-Tabs!)

R22	Harmful if swallowed.
R36	Irritating to eyes.
R36/38	Irritating to eyes and skin.
R41	Risk of serious damage to eyes.
R52/53	Harmful to aquatic organisms, may cause long term adverse effects in the aquatic environment.

### Further information:

This information is based on our present knowledge and describes the safety measures which are to be taken for this product. However, this shall not constitute a guarantee for any product features.

### Marking of changes:

A backslash in the left hand margin indicates an amendment from the previous version.