

TWO COLOR LINES

For the painted surfaces, you can choose neutral in white or add a stripe of semigloss metallic color to the water unit and the dentist's and assistant's elements.



THREE COLOR SCHEMES

When choosing your upholstery and paint colors, all options are open to you. To make your choice easier, we have grouped shades that go together particularly well in three color schemes.

THE ELEGANT WORLD
Classic, timeless combinations for an especially high-class look.

Upholstery colors:	Carbon	Aubergine	Platinum	Pacific
	Neutral white	Crimson red	Platinum silver	Polar blue



THE NATURAL WORLD
Soft, warm shades for a soothing treatment atmosphere.

Upholstery colors:	Basalt	Rosewood	Cherry	Sapphire
	Platinum silver	Metallic beige	Signal red	Polar blue



THE VITAL WORLD
Vivid colors and bright tones for a modern, lively atmosphere.

Upholstery colors:	Spring	Sunset	Orchid	Lagoon
	Fresh green	Dark orange	Pink	Turquoise blue



ALWAYS AT THE FOREFRONT OF INNOVATION!

As global innovation leader for dental equipment, we continuously invest in research and thus in the future of modern dentistry. By networking digital technologies with integrated solutions and optimizing the treatment workflow, we create improved treatment results, more comfort and safety for the patient as well as time and cost savings in everyday work. The combination of constant innovative power and globally growing sales and service structures makes Sirona the global market leader trusted by thousands of practices and labs around the world. **Enjoy every day. With Sirona.**



CAD/CAM systems
From pioneer to new standard. For 30 years we have been developing digital dentistry and creating new possibilities for the future practice and lab.

Imaging systems
Best image quality with the lowest dose. More than 100 years of developing x-rays for the dental practice make us the number 1 innovation partner.

Treatment centers
The business card of modern practices. We are striving to create the ideal ergonomic and innovative center. Individually tailored to the well-being and demands of the patient and dentist.

Instruments
Advantages that speak for themselves. We make sure that we provide the right balance of proven quality, individual ergonomics and innovative technology for user-friendly work.

Hygiene systems
Competence that gives you safety. When it comes to hygiene in the practice, we do not take any shortcuts.

YOU DEMAND EXCELLENCE. TENEO DELIVERS.



Subject to technical changes and errors in the text, Order No. A91100-1M41-8891-01-7600, Printed in Germany, Dispo-No. 04601, DEW15 WS 02/15.Y0

YOU WANT INNOVATION. TENEO BELONGS TO THE FUTURE.

The innovation class from Sirona is constantly evolving to continue to meet your needs tomorrow. This makes it a consistent choice for anyone who wants to work now and in the future with the most modern and comfortable technology while maintaining the highest standards of quality, design and workflow. With a variety of equipment options and the innovative, customizable control concept, TENEO is an investment in the future. **Enjoy every day. With Sirona.**



YOU WANT TO OFFER FIRST-CLASS TREATMENT. TENEO PUTS YOUR WORKFLOW IN MOTION.

Thanks to TENEO, you work at the highest level.

Integrated tabletop device functions, intuitive handling and innovative controls with many additional options allow you to offer treatment efficiently, quickly and comfortably. Benefit from incomparable comfort as a practitioner, while your patients are relaxed as never before.

HIGHEST POSSIBLE COMFORT FOR PROFESSIONAL AND PATIENT

Innovative workflow and automatic functions, ergonomic motor-driven sliding track, Thermo and Lounge Upholstery, lumbar support and massage functions

INTUITIVE USER INTERFACE

Easy to use, configurable and upgradable

INTEGRATION OF TABLE INSTRUMENT FUNCTIONS

For efficient endodontics and implantology treatment



YOU EXPECT THE BEST SUPPORT. TENEO IS RIGHT AT YOUR SIDE.

Ergonomically designed with smart integrated functions:
TENEO ensures healthy and comfortable work, so you can achieve perfect results from your first to your last patient.

Perfect headrest adjustment
The motorized headrest is automatically adjusted to the saved patient settings.*

Optimal working height
Regardless of the patient's stature: TENEO always adjusts to your suitable working height (vertical mouth point positioning).

NEW: Support arm system
High-quality materials and a functional design with uniform surfaces ensure long service life and easy cleaning.

Ergonomic working stools
The HUGO and CARL ergonomic working stools let you intuitively assume a correct, ergonomic working position.

NEW: Patient-specific chair positions
Patient-specific chair programs can be stored and recalled repeatedly.*

* Available with SIDEXIS 4.



Motor-driven sliding track
Programmable and operated via the foot control. This enables you to move the dentist element to the correct position effortlessly.




Optimal maneuverability
Flexible positioning of both dentist and assistant elements as well as independent tray positioning.



Innovative wireless foot control
For hands-free operation and free positioning with no cable clutter, allowing unrestricted movement sequences.



LEDview operating light
Together with the LED instrument light, optimum illumination of the entire treatment site is ensured.



YOU WANT INTUITIVE OPERATION. TENEO ADAPTS TO YOUR NEEDS.

The EasyTouch user interface is intuitive, individually programmable and enables the integration of new features at any time. With EasyTouch and the innovative Sirona instruments, you will deliver outstanding results even on long and busy working days.


NEW: Individually programmable quick-adjustment keys
Using quick-adjustment keys, you can set up your workflow to perform treatments in an optimal, time-saving manner.

Display
The display shows only those functions needed for the current treatment.


Instrumentation
Powerful motors, concentrated power, optimum speed, ergonomic weight distribution and perfect illumination. In this way, Sirona instruments ensure a relaxed treatment.

Hot keys
The hot keys are assigned to the functions that are regularly used.

With freely programmable settings and its selectable EasyMode, the EasyTouch menu guidance can very easily be adapted to your usual working methods. The endodontic and implantology functions can be integrated and thus make tabletop units and long setup times a thing of the past. And because the user interface is upgradable and the instruments are equipped with standardized interfaces, you can face the future with confidence.



Integrable implantology
With the integrable implant function, you can set the speed, torque and flow rate of the sterile cooling fluid individually to suit your workflow.



Selectable EasyMode
Shows all functions at a glance, ensuring easier, safer handling of the unit, and combines intuitive user-friendliness with the very highest treatment standards.

YOU WANT TO TREAT PATIENTS BOTH SAFELY AND EFFICIENTLY. TENEO GETS TO THE ROOT OF THINGS.

NEW:
EVEN MORE
ENDODONTICS
FUNCTIONS

Simplified workflow and perfect results: Thanks to the integrated and optimized endodontics function, your workflow is always optimal. You can work hygienically, safely, and efficiently.

The new and improved endodontic function guarantees a better workflow and thus more successful treatment. During treatment, you are aware of everything and can work without an additional tabletop unit or additional foot control.


Reciprocal file systems

Sirona is the only manufacturer of dental units to integrate the reciprocal file systems of Dentsply Maillefer and VDW directly into the operating system. Analogous to the tabletop units from these manufacturers, the reciprocal files can be selected and controlled as desired via the integrated file library from TENEO. Thanks to the reciprocal file systems, time-saving treatments with greater comfort and maximum safety are possible. Due to its one-time use, no disinfection, cleaning or sterilization is necessary, and cross-contamination is prevented.

wave•one®

wave•one®
GOLD

RECIPROC®
one file endo




Integrated endodontics

The integrated endodontics treatment function makes tabletop units a thing of the past and, along with precise torque control, also offers an optional ApexLocator for the highest degree of clinical safety. Other common file systems can be integrated in the unique file library.




Holder for the ApexLocator

The holder provides for hygienic storage of the ApexLocator during and after treatment. In addition, the practical positioning on the dentist element optimally supports your workflow.



ApexLocator on the monitor

When using the ApexLocator, the distance to the root canal end and the selected file can be displayed on the monitor. So you always have everything in sight.*



Contra-angle handpiece T1 LINE ENDO 6 L with light

Improved visibility during endodontic treatment increases the rate of successful treatment.

* Available with SIDEXIS 4.

YOU EXPECT SATISFIED PATIENTS. TENEO CREATES TRUST.

The innovation class stores patient-specific chair positions and provides your patients with maximum comfort in all respects. With the latest features to assist patient counseling, confidence is increased even further, ensuring long-term loyalty.

With TENEO, you have the latest tools for patient communication.. The practice network with Ethernet interface enables you to access intraoral images and other digital media directly from the console. You can use presentations to explain treatments, for example, or show video and x-ray images to illustrate treatment options to your patients.



Autofocus camera

The SiroCam AF* combines best image quality with easy operation. Easy automatic focusing and freeze frames by pressing the foot switch.



22" HD monitor

The HD monitor in 16:9 format provides best image quality – with integrated speakers as well as HDMI, DisplayPort and S-Video ports.



Lounge Upholstery

Soft, pleasant surface material and comfortable padding for utmost relaxation. Contoured surface for perfect positioning.

Thermo upholstery
Reduces the accumulated heat in the seat and back area and provides for more relaxation as a result of this cooling effect.

Massage function
The massage function in the backrest provides for additional relaxation during prolonged treatment.

Active lumbar support function
The back support of the anatomically shaped patient chair can be individually adapted to the spine of each individual patient.

Perfect positioning
Saved settings for each patient, unique chair elevation system, height adjustable from 370 to 820 mm for easy entry.



YOU WANT PERFECT RESULTS. TENEIO OFFERS MORE THAN ERGONOMICS.

Simplified workflow, perfect results:
Sirona is continually developing its ergonomics
program with this goal in mind.

Sirona's ergonomics program is based on decades of experience in ergonomics. And as an innovation driver for linking digital technologies with system solutions, we are constantly thinking further ahead. The ergonomics circle stands for the all-round vision of Sirona with regard to the connection between ergonomics and workflow in practice routine. Its round shape also indicates that Sirona's ergonomics program is being developed continuously.



INTUITIVE SITTING

More freedom of movement
Wireless foot control, motor-driven sliding track, high knee clearance under the backrest

Optimal maneuverability
Perfect coordination and flexible positioning of dentist and assistant elements, dentist's element and tray adjustable independently of each other

Intelligent working stools
HUGO, CARL

COMFORTABLE POSITIONING

Easy adjustment of the treatment center
Flexible height adjustment, optimal working height due to vertical mouth point positioning, patient-specific chair positions* and automatic headrest*

Best spine support
ErgoMotion, OrthoMotion, lumbar support function other

Greater patient comfort
Innovative chair elevation system, Premium thermo and Lounge Upholstery, massage function

OPTIMAL VISIBILITY

Best illumination
LEDview, LED instrument light

Motor-driven headrest
Automatically adjusts to current treatment step and patient*

Perfect intraoral images
Autofocus camera SiroCam AF*

INTEGRATED WORKFLOW

Intuitive operation
EasyTouch user interface with individualized operating concept

Integrated tabletop device functions
Implantology, endodontics with ApexLocator display on SIVISION monitor*

Modern equipment
Innovative instruments, modern patient communication system SIVISION, USB, and Ethernet interface

PERFECT RESULTS

Clinically perfect

Esthetically perfect

Perfect for the success of your practice



Click here for the ergonomics video at sirona.com

YOU WANT OPTIMAL HYGIENE. TENEO KEEPS YOU SAFE.

Hygiene in the practice has to be simple and absolutely safe.
With TENEO, you save time and still satisfy all requirements.

With TENEO, you can achieve the highest hygiene standards in minimum time. From the smooth, easy-to-clean surfaces to the easily detachable components and integrated sanitization adapter, everything is geared to the highest level of clinical safety. As a result, you can dispense with unnecessary work steps and take some of the workload off of your team.



Removable cuspidor

The motor-driven and thus hands-free cuspidor that can swivel in and out is easy to remove for easy cleaning and good practice hygiene.



Built-in disinfection unit with Aqua-Function

Permanent and automatic treatment of instrument spray and mouth rinsing water, sanitization of all critical water channels. Prevents manual operating errors. In case of an emergency, the Aqua-Function allows you to use the treatment center even without a connection to the water supply.



Thermo-disinfectable parts

All parts critical for hygiene are detachable and thermo-disinfectable or sterilizable. This increases protection from cross-contamination.



Built-in sanitization adapters

All instrument hoses can be connected to the integrated sanitization adapters for periodic cleaning of the waterlines. In addition, the suction hoses can be connected and cleaned from the inside with a stream of water at the press of a button, even during treatment. This function can optionally be done using chemicals from Sirona's list of cleaning and disinfectant agents.



The DGKH (German Society for Hospital Hygiene) has confirmed Sirona's adherence to the highest hygienic safety standards in accordance with the recommendations of the Robert Koch Institute.

YOU DEMAND QUALITY. TENEO SETS STANDARDS.

The starting point of the highest standards is the desire for absolute quality. The innovation class is therefore designed and manufactured with the highest quality standards in mind. With TENEO, you will work with long-term reliability and the highest level of quality.

Sophisticated ergonomics and high-quality materials help ensure an optimally designed workflow and create a foundation for perfect results. Whether you want to improve your workflow with the latest technology or integrate new developments: TENEO offers you greatest future sustainability.

Pioneering innovation

TENEO is continuously evolving so that you always work with the highest possible standard. The innovation class can be upgraded and, with a USB interface for integrating new functions and an Ethernet interface for easy troubleshooting by the service technician, it already has all the requirements today to assist you optimally in the future.

High-quality workmanship and design

For the innovation class, only the finest materials are used: Framework, chair and backrest support consist of a well-thought-out material combination of steel and aluminum. As a result, TENEO is light and yet extremely solid.

Quality Made in Germany

The innovation class is developed, manufactured and tested entirely at the Sirona factory in Bensheim. In our Center of Innovations, we are continually developing the product. As a result, the quality and reliability of your TENEO treatment unit are guaranteed over the long term – to ensure that your investment is protected.



YOU ARE LOOKING FOR THAT SOMETHING SPECIAL. **TENEO HAS STYLE.**

TENEO features a perfect combination of form and function. Prestigious awards such as the iF product design award and the nomination for the German Design award attest to the outstanding design of TENEO. With the three color schemes and our Lounge Upholstery, this treatment unit offers you a variety of design options for your practice. Be inspired and choose the configuration that is just right for you.



LOUNGE UPHOLSTERY.

The Lounge Upholstery not only affords a high degree of comfort, but also has an impressive dynamic design. The attractive contrast stitching and the continuous footrest give the Lounge Upholstery a classic look. The high-quality material has a soft, comfortable feel and is easy to clean and disinfect.



What will your treatment unit look like? Use our app to get TENE0 in your practice in advance and configure the desired upholstery color and furnishings right in your office.



YOUR INDIVIDUAL OPTIONS.






Options	Basic unit
Chair	
Patient-specific chair positions	■
Motor-driven headrest with automatic patient-specific headrest adjustment	■
Premium upholstery	■
Lounge Upholstery	□
Premium thermo upholstery	□
Lumbar support / massage function	□
Removable armrest, right armrest can be swiveled aside	□
Removable left armrest	□
4-way foot switch	■
Foot control, cable version	■
Wireless foot control	□
External device connection	□
Booster cushion for small patients	□
M1/C adapter plate, assembly without drilling	□
Dentist element (with EasyTouch, 6 positions, removable handles and instrument holder)	
SPRAYVIT M with light	□
1st drill drive (for motor and turbine)	■
2nd drill drive (for motor and turbine)	■
3rd drill drive (for motor and turbine)	■
4th drill drive (for motor and turbine)	□
BL electric motor	□
BL ISO C electric motor	□
BL implant electric motor	□
Turbine hose	■
SIROSONIC TL ultrasonic handpiece	□
SIROTOM electrosurgery kit	□
Implantology function	□
Endodontic treatment function	□
Reciprocating function	□
ApexLocator with display on SIVISION monitor	□
Integrated pump for saline solution	□
USB interface	■
Hands-free operation with cursor control	■
Contra-angle handpiece T1 LINE ENDO 6 L with light	□
Assistant's element (with 4 positions, removable holder)	
Satelec mini LED curing light	□
Surgical suction device	□
SPRAYVIT M with light	□

Options	Basic unit
Assistant element (with 4 positions, removable holder)	
Suction device	■
Saliva ejector	■
Suction hose with slide	□
Water unit	
Disinfection device with permanent disinfection and sanitization function	■
Suction hose cleaning with water	■
Suction hose cleaning with chemicals	□
Hydrocolloid port	□
Cuspidor valve for wet suction system	■
Automatic separator	□
Amalgam separator	□
Glass cuspidor, motor-driven swivel function, removable	■
Purge function with autopurge	■
Aqua-Function	■
SIVISION digital	
SiroCam AF*	□
22" monitor on the tray incl. tray	□
22" monitor on the light support column	□
Camera supply in the dentist element	□
Cable set for PC connection	□
Working stools (optional long or short gas spring)	
HUGO Freehand (with foot release)	□
HUGO Manual (with manual release)	□
HUGO/CARL Manual Plus (with manual release and foot ring)	□
CARL Manual (with manual release)	□
CARL Manual Plus (with manual release and foot ring)	□
Additional options	
Sirona Tray	□
PC interface (Ethernet)	■
PC control with EasyTouch	■
Remote diagnosis	■
LEDview operating light	□

* Available with SIDEXIS 4.
■ Basic unit or as a package. □ Optional.

YOUR EQUIPMENT OPTIONS.

DISCOUNT PACKAGES

COMFORT	INNOVATION	ENDO	IMPLANT	SIVISION
				
<ul style="list-style-type: none">2 HUGO or CARL working stools, selectable version and gas springThermo or optional Lounge UpholsteryLumbar support / massage function	<ul style="list-style-type: none">SPRAYVIT M for the dentist and assistant elements2 BL electric motors (upcharge for BL ISO C)SIROSONIC TL ultrasonic handpieceLEDview operating lightWireless foot control	<ul style="list-style-type: none">Endodontic treatment functionApexLocator functionT1 LINE ENDO 6 L contra-angle handpiece with lightISO adapter with light (optional)Holder for the ApexLocatorAdditional option: reciprocating function	<ul style="list-style-type: none">Implantology functionIntegrated pump for saline solutionBL implant electric motorIMPLANT 20:1 contra-angle handpiece	<ul style="list-style-type: none">22" monitorSiroCam AF*Monitor position on the tray or the light support columnCamera supply in the dentist's element



HUGO Manual Plus



HUGO Freehand



Booster cushion for small patients



Lounge Upholstery



NEW: Reciprocating function



External device connection



Integrated pump for saline solution



Hydrocolloid port



SiroCam AF*



NEW: HELIODENT^{PLUS} unit model



22" monitor on the tray



22" monitor on the light support arm