



SIROBoost – the new turbine with more power

The winners in the performance test are:
SIROBoost and You.

More power, more torque, more efficiency.

Power for your practice in the palm of your hand.

A dentist who does great things on a daily basis should also use powerful tools. The more powerful your instruments, the faster you can finish your procedures. And with SIROBoost, you can now have the most powerful turbine from Sirona. What does this power mean for you? Quite simply, you can work more efficiently right away with more torque.

Benefits at a glance:

- 22 watts of power for almost uninterrupted preparation
- Pleasant sound level, as SIROBoost is nice and quiet
- Makes your usual working routines even easier
- Optimal view with LED technology
- With even more sensation, you can save time and complete procedures faster
- You can concentrate on what is most important: your patient

First-class handheld technology.

The best thing for you about SIROBoost is: The powerful turbine. Time wasted stopping and starting work on a tooth is finally a thing of the past. With a powerful 22 watts, you now have even more torque and even more sensation. SIROBoost not only makes your usual routines easier, but also keeps the volume low. And the LED integrated in the coupling provides a perfect view of your field of treatment at all times. So you can concentrate fully on your patient.




High tech in high detail – that’s power.

The first LED to show its true colors.
Perfect vision with natural LED light.

Technical data	
Non-load speed (rpm)	approx. 420.000
Head diameter (mm)	11,4
Head height (mm)	14,8
Power (W)	approx. 22
Light/illuminance (lux)	approx. 25.000
Turbine weight (g)	approx. 57
Driving air pressure (bar)	2,7–3,0

- **Listen up!** Noise level significantly reduced for enhanced well-being
- **Let there be light!** 25,000 lux provide natural light intensity
- **That’s power!** A mighty 22 watts ensure maximum torque
- **The greatest of ease!** Ergonomic and anything but strenuous
- **Click&Go!** That’s how the fast change to Sirona works





Sirona – unique worldwide systems expertise in dental equipment products

Sirona develops and manufactures a comprehensive range of dental equipment, including CAD/CAM Systems for dental practices (CEREC) and laboratories (inLab), Instruments and Hygiene Systems, Treatment Centers and Imaging Systems. Sirona manufactures high technology products that guarantee ease of use and a high return on investment – for the good of your practice and for the benefit of your patients. In this way, you can approach every challenge you face with confidence. **Enjoy every day. With Sirona.**

Sirona Dental Systems · Fabrikstraße 31 · 64625 Bensheim
E-Mail: contact@sirona.com · www.sirona.com

The Dental Company

sirona.